For Binocular Microscope

LANOPTIK Technologies





www.lanoptik.com

Camera Solutions for Olympus Binocular Microscopes (1)



Functions Available:

	Solution 1	Solution 2	Solution 3
Deeply customized Android OS	\checkmark	\checkmark	
Pre-installed MS Office	\checkmark	\checkmark	
15.6-inch High Color Gamut Display	\checkmark		
Output Methods			
5G WiFi	\checkmark	\checkmark	\checkmark
USB			\checkmark
HDMI	\checkmark	\checkmark	\checkmark
WAN (Ethernet)			\checkmark

Camera Solutions for Olympus Binocular Microscopes (2)



Stereo Microscope

Applicable Model:

SZX10, etc.



Embedded Smart Camera



Dedicated for parallel optics type stereo microscope



Embedded WiFi Camera



Dedicated for parallel optics type stereo microscope

Functions Available:

	Solution 4	Solution 5
Deeply customized Android OS	\checkmark	
Pre-installed MS Office	\checkmark	
15.6-inch High Color Gamut Display		
Output Methods		
5G WiFi	\checkmark	\checkmark
USB		\checkmark
HDMI	\checkmark	\checkmark
WAN (Ethernet)		\checkmark

Camera Solutions for Nikon Binocular Microscopes (1)



I HDMI output is unavailable for 20MP cameras.

 \checkmark

USB HDMI

WAN (Ethernet)

Camera Solutions for Nikon Binocular Microscopes (2)



Stereo Microscope

Applicable Model:

SMZ 1270, SMZ800N, etc.



Embedded Smart Camera



Dedicated for parallel optics type stereo microscope



Embedded WiFi Camera



Dedicated for parallel optics type stereo microscope

Functions Available:

	Solution 4	Solution 5
Deeply customized Android OS	\checkmark	
Pre-installed MS Office	\checkmark	
15.6-inch High Color Gamut Display		
Output Methods		
5G WiFi	\checkmark	\checkmark
USB		\checkmark
HDMI	\checkmark	\checkmark
WAN (Ethernet)		\checkmark

Camera Solutions for Leica Binocular Microscopes (1)



Functions Available:

	Solution 1	Solution 2	Solution 3
Deeply customized Android OS	\checkmark	\checkmark	
Pre-installed MS Office	\checkmark	\checkmark	
15.6-inch High Color Gamut Display	\checkmark		
Output Methods			
5G WiFi	\checkmark	\checkmark	\checkmark
USB			\checkmark
HDMI	\checkmark	\checkmark	\checkmark
WAN (Ethernet)			\checkmark

Camera Solutions for Leica Binocular Microscopes (2)



Stereo Microscope

Applicable Model:

MZ10F, etc.



Dedicated for parallel optics type stereo microscope






Dedicated for parallel optics type stereo microscope

Functions Available:

	Solution 4	Solution 5
Deeply customized Android OS	\checkmark	
Pre-installed MS Office	\checkmark	
15.6-inch High Color Gamut Display		
Output Methods		
5G WiFi	\checkmark	\checkmark
USB		\checkmark
HDMI	\checkmark	\checkmark
WAN (Ethernet)		\checkmark

Camera Solutions for Zeiss Binocular Microscopes





LANOPTIK TECHNOLOGIES LTD

No. 72 Hongjing Street, Baiyun, Guangzhou, China. 510400Phone: +86 20 3898 6017 | Fax: +86 20 3847 6076Website: http://www.lanoptik.com | Email: info@lanoptik.com

